


Product Name: Recombinant Human AP4E1 Protein
Catalog #: PEH9524



Summary

Name	AP4E1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPM8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKEGYLPKKE SKTGDESGAL PVPQESIMEN VDQAITKKDQ SQVLTQSKEE KEKQLLASSL FVGLGSESTI

image



Background

Alternative Names	AP-4 complex subunit epsilon-1 (AP-4 adaptor complex subunit epsilon) (Adaptor-related protein complex 4 subunit epsilon-1) (Epsilon subunit of AP-4) (Epsilon-adaptin)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use