


Product Name: Recombinant Human APOL2 Protein
Catalog #: PEH9526



Summary

Name	APOL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9BQE5
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEGISFVLLD TGMGLGAAAA VAGITCSVVE LVNKLRRARQ ARNLDQSGTN VAKVMKEFVG GNTPNVLTIV DNWYQVTQGI GRNIRAIRRA RANPQLGAYA PPPHVIGRIS AEGGEQVERV VEGPAQAMSR GTMIVGAATG GILLLLDVVS LAYESKHLLE GAKSESAEEL KKRAQELEGK

image



Background

Alternative Names Apolipoprotein L2 (Apolipoprotein L-II) (ApoL-II)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human APOL2 Protein
Catalog #: PEH9526

