



Product Name: Recombinant Human YME1L1 Protein
Catalog #: PEH9529



Summary

Name	YME1L1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96TA2
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTGLDSAVDP VQMKNVTFEH VKGVVEAKQE LQEVVEFLKN PQKFTILGGK LPKGILLVGP PGTGKTLLAR AVAGEADVFP YYASGSEFDE MFVGVGASRI RNLFREAKAN APCVIFIDEL DSVGGKRIES PMHPYSRQTI NQLLAEMDGF KPNEGVIIG ATNFPEALDN ALIRPGRFDM QVTVPRPDVK GRTEILKWYL NLIKFDQSVD PEIARGTVG FSGAELENLV NQAALKAARD GKEMVTMKEL EFSKDILMG PERRSVEIDN KNKTITAYHE SGHAIAYYT KDAMPINKAT IMPRGPTLGH VSLLPENDRW NETRAQLLAQ MDVSMGGRVA EELIFGTDHI TTGASSDFDN ATKIAKRMVT KFGMSEKLGV MTYSDTGKLS PETQSAIEQE IRILLRDSYE RAKHILKTHA KEHKNLAEAL LTYETLDAKE IQIVLEGKKL

image



Background

Product Name: Recombinant Human YME1L1 Protein
Catalog #: PEH9529



Alternative Names

ATP-dependent zinc metalloprotease YME1L1 (EC 3.4.24.-) (ATP-dependent metalloprotease FtsH1) (Meg-4) (Presenilin-associated metalloprotease) (PAMP) (YME1-like protein 1)

Note

For Research Use Only , Not for Diagnostic Use