


Product Name: Recombinant Human BCO1 Protein
Catalog #: PEH9534



Summary

Name	BCDO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAY6
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DGVILSAIVS TDPQKLPELL ILDAKSFTSL ARASVDVDMH MDLHGLFITD MDWDTKKQAA SEEQRDRASD

image



Background

Alternative Names	Beta,beta-carotene 15,15'-dioxygenase (EC 1.13.11.63) (Beta-carotene dioxygenase 1) (Beta-carotene oxygenase 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use