


Product Name: Recombinant Human KCTD9 Protein
Catalog #: PEH9541



Summary

| | |
|------------------------------|--|
| Name | KCTD9_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q7L273 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 43 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SPKNGKVVAV YGTLSDLLSV ASSKLGIKAT SVYNGKGGLI DDIALIRDDD VLFVCEGEPF IDPQTDSKPP ELLGFHTDW LTLNVGGRYF TTRSTLVNK EPDSMLAHMF KDKGVWGNKQ DHRGAFLIDR SPEYFEPILN YLRHGQLIVN DGINLLGVLE EARFFGIDSL IEHLEVAIKN SQPPEDHSPI SRKEFVRFL ATPTKSELRC QGLNFSGADL SRLDLRYINF KMANLSRCNL AHANLCCANL ERADLSGSVL DCANLQGVKM LCSNAEGASL KLCNFEDPSG LKANLEGANL KGVDMEGSQM TGINLRVATL KNAKLKNCNL RGATLAGTDL ENCDLSGCDL QEANLRGSNV KGAIFFEMLT |

image



Background

Alternative Names BTB/POZ domain-containing protein KCTD9

Product Name: Recombinant Human KCTD9 Protein
Catalog #: PEH9541



Note

For Research Use Only , Not for Diagnostic Use