



Summary

Name	CAYP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXY5
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FSYIPLSTAN LPNSSSTLGW VTPCQTPYTQ YHLNKLDQNI IPENLPAPTD KCKLKYQQCK TEIKEGYKQY SQRNAENTKS NVTHKQSPRN KIDEKCVQDE EANTDDLTTL DRKAILQQGY ADNSCDKQQR ARKLD AEIVA AEKKKQIVAE QVMIDHLSRA VISDPEQNLA IEQKESDHIL PDSKMTPLRF RKRTLHETKI RTHSTLTENV LSHKLQFDGR

image



Background

Alternative Names Calcyphosin-2 (Calcyphosine-2)

Note

Product Name: Recombinant Human CAPS2 Protein
Catalog #: PEH9542



For Research Use Only , Not for Diagnostic Use