


**Product Name: Recombinant Human CAPN10 Protein**  
**Catalog #: PEH9543**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | CAN10_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | N-SUMO & C-His  |
| <b>Accession #</b>           | Q9HC96  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 48 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | RELFRDAAFP AADSSLFCDL STPLAQFRED ITWRRPQEIC ATPRLFDDP<br>REGQVKQGLL GDCWFLCACA ALQKSRHLLD QVIPPGQPSW ADQEYRGSFT<br>CRIWQFGRWV EVTTDDRLPC LAGRLCFSRC QREDFVWLPL LEKVYAKVHG<br>SYEHLWAGQV ADALVDLTGG LAERWNLKGV AGSGGQQDRP GRWEHRTCRC<br>LLHLKDQCLI SCCVLSPRAG ARELGEFHAF IVSDLRELQG QAGQCILLR<br>IQNPWGRRCW QGLWREGGEG WSQVDAAVAS ELLSQLQEGE FWVEEEFLR<br>EFDELTVGYF |

## image



## Background

**Alternative Names** Calpain-10 (EC 3.4.22.-) (Calcium-activated neutral proteinase 10) (CANP 10)

**Product Name: Recombinant Human CAPN10 Protein**  
**Catalog #: PEH9543**

---



**Note**

For Research Use Only , Not for Diagnostic Use