


Product Name: Recombinant Human CT55 Protein
Catalog #: PEH9546



Summary

Name	CT55_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WUE5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQGPQQQLP QGDTQLTTVQ GVVTSFCGDY GMIDESIYFS SDVVTGNVPL KVGQKVNVVV EEDKPHYGLR AIKVDVPRH LYGAGPSDSG

image



Background

Alternative Names Cancer/testis antigen 55 (Tumor antigen BJ-HCC-20)

Note

For Research Use Only , Not for Diagnostic Use