


Product Name: Recombinant Human CEP55 Protein
Catalog #: PEH9552



Summary

Name	CEP55_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q53EZ4
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSKWGSKPSN SKSETTLEKL KGEIAHLKTS VDEITSGK GK LTDKERHRL EKIRVLEAEK EKNAYQLTEK DKEIQRLRDQ LKARYSTTTL LEQLEETTRE GERREQVLKA LSEEKDVLKQ QLSAATSRIA ELESKTNTLR LSQTVAPNCF NSSINNIHEM EIQLKDALEK NQQWLVYDQQ REVYVKGLLA KIFELEKKTE TAAHSLPQQT KKPESEGYLQ EEKQKCYNDL LASAKKDLEV ERQTITQLSF ELSEFRRKYE ETQKEVHNLN QLLYSQRRAD VQHLEDDRHK TEKIQLKREE NDIARGKLEE

image



Background

Alternative Names Centrosomal protein of 55 kDa (Cep55) (Up-regulated in colon cancer 6)

Product Name: Recombinant Human CEP55 Protein
Catalog #: PEH9552



Note

For Research Use Only , Not for Diagnostic Use