


Product Name: Recombinant Human CEP83 Protein
Catalog #: PEH9553



Summary

Name	CEP83_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y592
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RILDEGKIKY ESEIARLEED KEELRNQLLN VDLTKDSKRV EQLAREKVYL CQKLKGLEAE VAEKKAEN SEAQVENAQR IQVRQLAEMQ ATVRSLEAEK QSANLRAERL EKELQSSSEQ NTFLINKLHK AEREINTLSS KVKEKHSNK LEITDIKLET ARAKSELERE RNKIQSELDG LQSDNEILKA AVEHHKVLLV EKDRELIRKV QAAKEEGYQK LVVLQDEKLE LENRLADLEK

image



Background

Alternative Names	Centrosomal protein of 83 kDa (Cep83) (Coiled-coil domain-containing protein 41) (Renal carcinoma antigen NY-REN-58)
--------------------------	--

Note

Product Name: Recombinant Human CEP83 Protein
Catalog #: PEH9553



For Research Use Only , Not for Diagnostic Use