

Product Name: Recombinant Human CERK Protein
Catalog #: PEH9554



Summary

Name	CERK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TCT0
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EMLEKLTSRP KHLVFINPF GGKGGKRIY ERKVAPLFTL ASITTDIIVT EHANQAKETL YEINIDKYDG IVCVGGDGMF SEVLHGLIGR TQRSAGVDQN HPRAVLVPSS LRIGIIPAGS TDCVCYSTVG TSDAETSALH IVVGDSLAMD VSSVHHNSTL LRYSVSLLGY GFYGDIIKDS

image



Background

Alternative Names Ceramide kinase (hCERK) (EC 2.7.1.138) (Acylsphingosine kinase) (Lipid kinase 4) (LK4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CERK Protein
Catalog #: PEH9554

