


Product Name: Recombinant Human CDR2L Protein
Catalog #: PEH9555



Summary

Name	CDR2L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86X02
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAEEEEWYD QQDLEQDLHL AAELGKTLLE RNKELEGLSQ QMYSTNEEQV QEIEYLTQQL DTLRHVNEQH AKVYEQLDLT ARDLELTNHR LVLESKAAQQ KIHGLTETIE RLQAQVEELQ AQVEQLRGLQ QLRVLREKRE

image



Background

Alternative Names	Cerebellar degeneration-related protein 2-like (Paraneoplastic 62 kDa antigen)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use