



Product Name: Recombinant Human CFAP45 Protein
Catalog #: PEH9557



Summary

Name	CFA45_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UL16
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HTLQKTLTAL GLDRKPETIQ LITRDMVREL IVPTEDPSGE SLIISPEEFE RIKWASHVLT REELEARDQA FKKEKEATMD AVMTRKKIMK QKEMVWNNNK KLSDLEEVAK ERAQNLLQRA NKL RMEQEEE LKDMSKIILN AKCHAIRDAQ ILEKQQIQKE LDTEEKRLDQ MMEVERQKSI QRQEELERKR REERIRGRRQ IVEQMEKNQE ERSLLAEQRE QEKEQMLEYM EQLQEEDLKD MERRQQQKLK MQAEIKRIND ENQKQKAELL AQEKLADQMV MEFTKKKMAR EAEFEAEQER IRREKEKEIA RLRAMQEKAQ DYQAEQDALR AKRNQEVADR EWRRKEKENA RKKMETEAE L RKSRLQVAF KEHALAVQVQ RDRDEFERIL RAQREQIEKE

image



Background

Alternative Names	Cilia- and flagella-associated protein 45 (Coiled-coil domain-containing
--------------------------	--

Product Name: Recombinant Human CFAP45 Protein
Catalog #: PEH9557



protein 19) (Nasopharyngeal epithelium-specific protein 1)

Note

For Research Use Only , Not for Diagnostic Use