


Product Name: Recombinant Human CLHC1 Protein
Catalog #: PEH9560



Summary

Name	CLHC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NHS4
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEQIQKTKYL QGDQGIVEEL MEDDPRRAKE AEIMLHYIER FNELISLGEY EKAACYAANS PRRILRNIGT MNTFKAVGKI RGKPLPLLLF FEALFITSHA FPCPVDAALT LEGIKCGLSE KRLDLVTNWV TQERLTFSEE AGDVICDYGE QDTYNKAKCL ALAQIVYSEC GLHKKAILCL CKQGQTHRVM EYIQQLKDFT TDDLQLLMS CPQVELIQCL TKELNEKQPS LSFGLAILHL FSADMKKVGI KLLQEINKGG IDAVESLMIN DSFCSIEKWQ EVANICSQNG FDKLSNDITS ILRSQAAVTE

image



Background

Alternative Names Clathrin heavy chain linker domain-containing protein 1

Product Name: Recombinant Human CLHC1 Protein
Catalog #: PEH9560



Note

For Research Use Only , Not for Diagnostic Use