


Product Name: Recombinant Human COPG2 Protein
Catalog #: PEH9561



Summary

Name	COPG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBF2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEITLVATKP EKLAPSRQDI FQEQLAAIPE FLNIGPLFKS SEPVQLTEAE TEYFVRCIKH MFTNHIVFQF DCTNTLNDQL LEKVTVMQMEP SDSYEVLSCL PAPSLPYNQP GICYTLVRLP DDDPTAVAGS FSCTMKFTVR DCDPNTGVDP EDGYDDEYVL EDLEVTVDH IQKVLKPNFA AAWEEVGDTF EKEETFALSS TKTLEEAVNN IITFLGMQPC

image



Background

Alternative Names	Coatmer subunit gamma-2 (Gamma-2-coat protein) (Gamma-2-COP)
--------------------------	--

Note

Product Name: Recombinant Human COPG2 Protein
Catalog #: PEH9561



For Research Use Only , Not for Diagnostic Use