


Product Name: Recombinant Human CCDC170 Protein
Catalog #: PEH9564



Summary

Name	CC170_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IYT3
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VELNEKLQKC SKENEENKKQ VSKNCRKHEE FLTQLRDCLD PDERNDKASD EDLILKLRDL RKENEFVKGQ IVILEETINV HEMEAKASRE TIMRLASEVN REQKKAASCT EEKEKLNQDL LSAVEAKEAL EREVKIFQER LLAGQQVWDA SKQEVSLKK SSSELEKSLK ASQDAVTTSQ SQYFSFREKI AALLRGRISM TGSTEDTILE KIREMDSREE SRDRMVSQL E AQISELVEQL GKESGFHQKA LQRAQKAENM LETLQGQLTH LEAELVSGGV LRDNLNFEKQ KYLKFLDQLS QKMKLDQMAA ELGFDMLRDV VLARTEQLVR LESNAVIENK TIAHNLQRKL KTQKERLESK ELHMSLLRQK IAQLEEEKQA RTALVVERDN AHLTIRNLQK KVERLQKELN TCRDLHTELK AKLADTNELK IKTLEQTKAI EDLNKSRDQL EKMKEKAEEK LMSVKSELD TTEHEAKENKE

image



Background

Product Name: Recombinant Human CCDC170 Protein
Catalog #: PEH9564



Alternative Names

Coiled-coil domain-containing protein 170

Note

For Research Use Only , Not for Diagnostic Use