


Product Name: Recombinant Human CCDC28A Protein
Catalog #: PEH9565



Summary

Name	CC28A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWP9
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVNAKKNAIP VSKSTGFSNP ASQSTSQRPK LKRVMKEKTK PQGGEKGGAQ STPIQHSFLT DVSDVQEMER GLLSLLNDFH SGKLAQFGNE CSIEQMEHVR GMQEKLARLN LELYGELEEL PEDKRKTASD SNLDRLLSDL EELNSSIQKL HLADAQDVPN

image



Background

Alternative Names Coiled-coil domain-containing protein 28A (CCRL1AP)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CCDC28A Protein
Catalog #: PEH9565

