


Product Name: Recombinant Human CCDC30 Protein
Catalog #: PEH9566



Summary

Name	CCD30_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VVM6
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELETKVLKLS QEFAQLNHFT LGGKTAPSNL ITSENTCKDP ESNEPILETE IQSRKEETEE LCPKLGERKQ KEIPEESVKE GSFPRREGQKE EGSQQNRDMK DEEKEQQLTM KPVEIVRLRE ELSHINQSL QSQSSGDSSD DSGAQHPSSG EKLKYNQQGE VQQLHQNLHR LQILCNSAEN ELRYERGQNL DLKQHNSLLQ EENIKIKIEL KHAQQKLLDS TKMCSSLTAE YKHCQQKIKE LELEVLKHTQ SIKSQNNLQE KLVQEKSVA DAEKILDLQ RKLEHAHKVC LTDTCISEKQ QLEEKIKEAT QNEAKVKQQY QEEQQKRKLL YQNVDELHRQ VRTLQDKENL LEMTCSQQQS RIQQQEALLK QLENEKRKYD EHVKSQNELS EKLSKLQKEK EALREEYLRL LKLLNVHVRN YNEKHHQQKV KLQVKYRLT NEVELRDKRI NQFEDEIGIL QHKIEKEKAI QDQITAQNDT

image



Background

Product Name: Recombinant Human CCDC30 Protein
Catalog #: PEH9566



Alternative Names

Coiled-coil domain-containing protein 30 (Prefoldin subunit 6-like protein)

Note

For Research Use Only , Not for Diagnostic Use