


Product Name: Recombinant Human CCDC83 Protein
Catalog #: PEH9569



Summary

Name	CCD83_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWF9
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTHDGPPKEI KLPTSEALLD YQCQIKEDAV EQFMFQIKTL RKKNQKYHER NSRLKEEQIW HIRHLLKELS EEKAEGLPVV TREDVEEAMK EKWKFERDQE KNLRDMRMQI SNAEKLFEK LSEKEYWEEY KNVGSRHAK LITSLQNDIN TVKENAEKMS EHYKITLEDT RKKIHKETLL QLDQKKEWAT QNAVKLIDKG SYLEIWENDW LKKEVAIHRK EVEELKNAIH ELEAENLVLI DQLSNCRLVD LKIPRRLYLT QAAGLEVPPE EMSLELPETH IEEKSELQPT EVESRDLMSS SDESTILHLS HENSIEDLQY VKIDKEENSG TFGDITDMKY LLYEDEKDFK DYVNLGPLGV KLMSVESKKM PIHFQEKEIP

image



Background

Alternative Names	Coiled-coil domain-containing protein 83
--------------------------	--

Product Name: Recombinant Human CCDC83 Protein
Catalog #: PEH9569



Note

For Research Use Only , Not for Diagnostic Use