


Product Name: Recombinant Human CCDC93 Protein
Catalog #: PEH9570



Summary

Name	CCD93_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q567U6
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KYKRHQGAEE LLDEESRIHA TLLEYGRRYG FSRQSKMEKA EDKKTALPAG LSATEKADAH EEDELRAAEE QRIQSLMTKM TAMANEESSL TASSVGQIVG LCSAEIKQIV SEYAEKQSEL SAEESPEKLG TSQLHRRKVI SLNKQIAQKT KHLEELRASH TSLQARYNEA KKTLELTKY SEKLDKEQAA LEKIESKADP SILQNLRALV AMNENLKSQE QEFKAHCREE MTRLQQEIEN LKAERAPRGD EKTLSGGEPP GTLTSAMTHD EDLDRRYNME

image



Background

Alternative Names Coiled-coil domain-containing protein 93

Product Name: Recombinant Human CCDC93 Protein
Catalog #: PEH9570



Note

For Research Use Only , Not for Diagnostic Use