


Product Name: Recombinant Human CCDC96 Protein
Catalog #: PEH9571



Summary

Name	CCD96_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2M329
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETERVEGEEE DKEETQRDGA ESKERDGEGR PAKSQEEGKR LYGRDEFEDL EWSEEVQKLQ EQQLRSDLLD QYRSLVERN RSQRYNLYLQ HKIFEALRRK KGLEAAEVAD RGAEAEAPEK EQAYLRHLGM LEELKKQQAD DLQWYHQELG QLKRQCQEKL TRVEKEWRRF QALKKQVVMQ AMGSCRMGG RQAALREVEQ IQALEDKKEK EMSAVRENI QLKQSLVHFE

image



Background

Alternative Names Coiled-coil domain-containing protein 96

Note

Product Name: Recombinant Human CCDC96 Protein
Catalog #: PEH9571



For Research Use Only , Not for Diagnostic Use