


Product Name: Recombinant Human NCAPG Protein
Catalog #: PEH9574



Summary

Name	CND3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BPX3
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RTCLRALEKI KIQLEKGNKE FGDQAEAAQD ATLTTTTFQN EDEKNKEVYM TPLRGVKATQ ASKSTQLKTN RGQRKVTVSA RTNRRCQTAE ADSESDHEVP EPESEMKMRL PRRAKTAALE

image



Background

Alternative Names	Condensin complex subunit 3 (Chromosome-associated protein G) (Condensin subunit CAP-G) (hCAP-G) (Melanoma antigen NY-MEL-3) (Non-SMC condensin I complex subunit G) (XCAP-G homolog)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NCAPG Protein
Catalog #: PEH9574

