


Product Name: Recombinant Human CPNE5 Protein
Catalog #: PEH9575



Summary

Name	CPNE5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCH3
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FNIYEVVNP KKMKKKKYVN SGTVTLLSFA VESECTFLDY IKGQTQINFT VAIDFTASNG NPSQSTSLHY MSPYQLNAYA LALTAVGEII QHYDSDKMFP ALGFGAKLPP DGRVSHEFPL NGNQENPSCC GIDGILEAYH RSLRTVQLYG PTNFAPVVTH VARNAAAVQD GSQYSVLLII TDGVISDMAQ TKEAIVNAAK LPMSIIIVGV GQAEFDAMVE LDGDDVRISS RGKLAERDIV QFVPFRDYVD RTGNHVL SMA RLARDVLA EI PDQLVSY MKA QGIRPRPPPA

image



Background

Alternative Names Copine-5 (Copine V)

Product Name: Recombinant Human CPNE5 Protein
Catalog #: PEH9575



Note

For Research Use Only , Not for Diagnostic Use