


Product Name: Recombinant Human MMAB Protein
Catalog #: PEH9576



Summary

Name	MMAB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EY8
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDGDRPQSS KTRIPKIYT KTGDKGFSST FTGERRPKDD QVFEAVGTTD ELSSAIGFAL ELVTEKGHTF AEELQKIQCT LQDVGSALAT PCSSAREAHL KYTTFKAGPI LELEQWIDKY TSQLPPLTAF ILPSGGKISS ALHFCRAVCR RAERRVVPLV QMGETDANVA KFLNRLSDYL FTLARYAAMK EGNQEKIYMK

image



Background

Alternative Names	Corrinoid adenosyltransferase (EC 2.5.1.17) (Cob(II)alamin adenosyltransferase) (Cob(II)yrinic acid a,c-diamide adenosyltransferase) (Cobinamide/cobalamin adenosyltransferase) (Methylmalonic aciduria type B protein)
--------------------------	---

Product Name: Recombinant Human MMAB Protein
Catalog #: PEH9576



Note

For Research Use Only , Not for Diagnostic Use