


Product Name: Recombinant Human CELF4 Protein
Catalog #: PEH9577



Summary

Name	CELF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9BZC1
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GNPSTIPMKD HDAIKLFIGQ IPRNLDEKDL KPLFEEFGKI YELTVLKDRF TGMHKGCAFL TYCERESALK AQSALHEQKT LPGMNRPIQV KPADSESRRG SSCLRQPPSQ

image



Background

Alternative Names	CUGBP Elav-like family member 4 (CELF-4) (Bruno-like protein 4) (CUG-BP-and ETR-3-like factor 4) (RNA-binding protein BRUNOL-4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use