

Product Name: Recombinant Human CXXC4 Protein
Catalog #: PEH9579



Summary

Name	CXXC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2H0
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAFQIANLAD CPQNHSSSS SSSGGAGGAN PAKKKRKRCG VCVPCRRLIN CGVCSSCRNR KTGHQICKFR KCEELKKKPG TSLERTPVPS

image



Background

Alternative Names	CXXC-type zinc finger protein 4 (Inhibition of the Dvl and axin complex protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use