


Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9582



Summary

Name	CSAD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9Y600
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEAEGAVPFL VSATSGTTVL GAFDPLEAIA DVCQRHGLWL HVDAAWGGSV LLSQTHRHL DGIQRADSV WNPBKLLAAG LQCSALLQD TSNLLKRCHG SQASYLFQQD KFYDVALDTG DKVVQCGRRV DCLKLWLMWK AQGDQGLERR IDQAFVLARY LVEEMKKREG

image



Background

Alternative Names	Cysteine sulfinic acid decarboxylase (EC 4.1.1.29) (Aspartate 1-decarboxylase) (EC 4.1.1.11) (Cysteine-sulfinic acid decarboxylase) (Sulfinic acid decarboxylase)
--------------------------	---

Note

Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9582



For Research Use Only , Not for Diagnostic Use