


Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9584



Summary

Name	CSAD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9Y600
Host	E.coli
Species	Human
Molecular Mass	87 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGDPVAVEAL LRAVFGVVVD EAIQKGTSVS QKVCEWKEPE ELKQLLDLEL RSQGESSQKQI LERCRAVIRY SVKTGHPREF NQLFSGLDPH ALAGRIITES LNTSQYTYEI APVFVLMEEE VLRKLRLALVG WSSGDGIFCP GGSISNMYAV NLARYQRYPD CKQRGLRTLPLALFTSKEC HYSIQKGAAP LGLGTDSVRV VKADERGKMV PEDLERQIGM AEAEGAVPFL VSATSGTTVL GAFDPLEAIA DVCQRHGLWL HVDAAWGGSV LLSQTHRHL DGIQRADSV WNPBKLLAAG LQCSALLQD TSNLLKRCHG SQASYLFQQD KFYDVALDTG DKVVQCGRRV DCLKLWLMWK AQGDQGLERR IDQAFVLARY LVEEMKKREG

image



Background

Alternative Names Cysteine sulfinic acid decarboxylase (EC 4.1.1.29) (Aspartate 1-

Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9584



decarboxylase) (EC 4.1.1.11) (Cysteine-sulfinate decarboxylase)
(Sulfinioalanine decarboxylase)

Note

For Research Use Only , Not for Diagnostic Use