


Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9585



Summary

Name	CSAD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9Y600
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQASYLFQQD KFYDVALDTG DKVVQCGRRV DCLKLWLMWK AQGDQGLERR IDQAFVLARY LVEEMKKREG FELVMEPEFV NVCFWFVPPS LRGKQESPDY HERLSKVAPV LKERMVKEGS MMIGYQPHGT RGNFFRVVVA NSALTCADMD FLLNELERLG

image



Background

Alternative Names	Cysteine sulfinic acid decarboxylase (EC 4.1.1.29) (Aspartate 1-decarboxylase) (EC 4.1.1.11) (Cysteine-sulfinic acid decarboxylase) (Sulfinic acid decarboxylase)
--------------------------	---

Note

Product Name: Recombinant Human CSAD Protein
Catalog #: PEH9585



For Research Use Only , Not for Diagnostic Use