


Product Name: Recombinant Human CIAO2A Protein
Catalog #: PEH9587



Summary

Name	CIA2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H5X1
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LSEPGAARQP RIMEEKALEV YDLIRTIRD P EKPNTLEELE VVSESCVEVQ EINEEEYLVI IRFTPTVPHC SLATLIGLCL RVKLQRCLPF KHKLEIYISE GTHSTEEDIN KQINDKERVA AAMENPNLRE IVEQCVLEPD

image



Background

Alternative Names Cytosolic iron-sulfur assembly component 2A (MIP18 family protein FAM96A)

Note

For Research Use Only , Not for Diagnostic Use