


Product Name: Recombinant Human DIP2B Protein
Catalog #: PEH9593



Summary

Name	DIP2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P265
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAVAALPPEV RAQLAELELE LSEGDIQKG YEKKRSKLLS PYSPQTQETD SAVQKELRNQ TPAPSAAQTS APSKYHRTRS GGARDERYRS DIHTEAVQAA LAKHKEQKMA LPMPKRRST FVQSPADACT PPDTSASED EGSLRRQAAL SAALQQSLQN AESWINRSIQ GSSTSSSASS TLSHGEVKGT SGLADVFAN TRIEFNSAPP DVTTTTSSSS SSSSIRPANI DLPPSGIVKG MHKGSNRSSL

image



Background

Alternative Names Disco-interacting protein 2 homolog B (DIP2 homolog B)

Note

Product Name: Recombinant Human DIP2B Protein
Catalog #: PEH9593



For Research Use Only , Not for Diagnostic Use