


Product Name: Recombinant Human DUSP13 Protein
Catalog #: PEH9600



Summary

Name	DS13B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UII6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDSLQKQDLR RPKIHGAVQA SPYQPPTLAS LQRLWVRQA ATLNHIDEVW PSLFLGDAYA ARDKSKLIQL GITHVVNAAA GKFQVDTGAK FYRGMSLEY GIEADDNPFF DLSVYFLPVA RYIRAALSVP QGRVLVHCAM GVSRSATLVL AFLMICENMT LVEAIQTVQA HRNICPNSGF LRQLQVLDNR LGRETGRF

image



Background

Alternative Names	Dual specificity protein phosphatase 13 isoform B (DUSP13B) (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase SKRP4) (Testis- and skeletal-muscle-specific DSP)
--------------------------	---

Note

Product Name: Recombinant Human DUSP13 Protein
Catalog #: PEH9600



For Research Use Only , Not for Diagnostic Use