


Product Name: Recombinant Human CDC14A Protein
Catalog #: PEH9602



Summary

Name	CC14A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNH5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QHFLEEKQAS LWVQGDIFRS KLKNRPSSEG SINKILSGLD DMSIGGNLSK TQNMERFGED NLEDDDVEMK NGITQGDKLR

image



Background

Alternative Names	Dual specificity protein phosphatase CDC14A (EC 3.1.3.16) (EC 3.1.3.48) (CDC14 cell division cycle 14 homolog A)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use