


Product Name: Recombinant Human DNAI3 Protein
Catalog #: PEH9603



Summary

Name	DNAI3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWG1
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGKKRLKPVL AASEDMEPVN MESMGHPEIY PLVLTTKTQE IFNCRIDEDV TDEQPYKLIN KEDIFEDLRN RAAVSDFHPV KKIVQEYPGN ELLVYDKDF KYGLNFYLIA TEEGKENYLN PPEVP EEQEE YKEHIPEDVY IYKPPVSKPW VSLGSEKEIE EESVTESTKQ ITYMISRKRS EFGAPIKFSD QNASSVKDAY IECTAYPDKN FTLKQLEKDV GMQVIPQIKD ISTQTKWTYP KNATTQYYPR EFSEEEKETL KQSKPLVDFL NNASISVEIA LQQNEIMNTF IDDWKYLAEE EGTFGDKTDT HLKEYQSFTD LHSPTKMIT

image



Background

Alternative Names	Dynein intermediate chain 3, axonemal (Dynein axonemal intermediate chain 3) (Testis development protein NYD-SP29) (WD repeat-containing
--------------------------	--

Product Name: Recombinant Human DNAI3 Protein
Catalog #: PEH9603



protein 63)

Note

For Research Use Only , Not for Diagnostic Use