


Product Name: Recombinant Human MARCHF7 Protein
Catalog #: PEH9608



Summary

Name	MARH7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H992
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLQKIKESLL LEDSEEEEGD LCRICQMAAA SSSNLLIEPC KCTGSLQYVH QDCMKKWLQA KINS GS SLEA VTTCELCKEK LELNLEDFDI HELHRAHANE QAEYEFISSG LYLVLHLHC EQSFSDMMGN TNEPSTRVRF INLARTLQAH MEDLETSEDD

image



Background

Alternative Names	E3 ubiquitin-protein ligase MARCHF7 (EC 2.3.2.27) (Axotrophin) (Membrane-associated RING finger protein 7) (Membrane-associated RING-CH protein VII) (MARCH-VII) (RING finger protein 177) (RING-type E3 ubiquitin transferase MARCHF7)
--------------------------	---



Product Name: Recombinant Human MARCHF7 Protein
Catalog #: PEH9608

Note

For Research Use Only , Not for Diagnostic Use