


Product Name: Recombinant Human MSL2 Protein
Catalog #: PEH9609



Summary

Name	MSL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCI7
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVQPLISTI IRGPTLGASA PVTVKRESKI SLQPIATVPN GGTPKISKT VLLSTKSMKK SHEHGSKKSH SKTKPGILKK DKAVKEKIPS HHFMPGSPTK TVYKKPQEEK

image



Background

Alternative Names	E3 ubiquitin-protein ligase MSL2 (EC 2.3.2.-) (E3 ubiquitin-protein transferase MSL2) (Male-specific lethal 2-like 1) (MSL2-like 1) (Male-specific lethal-2 homolog) (MSL-2) (Male-specific lethal-2 homolog 1) (RING finger protein 184)
--------------------------	---

Note

Product Name: Recombinant Human MSL2 Protein
Catalog #: PEH9609



For Research Use Only , Not for Diagnostic Use