


Product Name: Recombinant Human PELI3 Protein
Catalog #: PEH9610



Summary

Name	PELI3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N2H9
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPRTSDLQHR GNKGSCVLSS PGEDAQPGEE PIKYGELIVL GCCEEGGEET EAQRGEVTGP RAHSCYNGCL ASGDKGRRRS RLALSRRSHA NGVKPDVMHH ISTPLVSKAL SNRGQHSISY TLRSHSVIV EYTHDSDTDM FQIGRSTENM IDFVVTDTSPP GGGAAEGPSA

image



Background

Alternative Names	E3 ubiquitin-protein ligase pellino homolog 3 (Pellino-3) (EC 2.3.2.27) (RING-type E3 ubiquitin transferase pellino homolog 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PELI3 Protein
Catalog #: PEH9610

