

Product Name: Recombinant Human SH3RF3 Protein
Catalog #: PEH9612



Summary

Name	SH3R3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEJ3
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AHSAAQAQDR PTATVSPLRT QNSPSRLPAT SLRPHSVVSP QHSHQPPVQM CPRPAIPLTS AASAITPPNV SAANLNGEAG GGPIGVLSTS SPTNTGCKLD EKKSEKKEKK SGLLKLLAGA STKKKSRSP

image



Background

Alternative Names	E3 ubiquitin-protein ligase SH3RF3 (EC 2.3.2.27) (Plenty of SH3s 2) (SH3 domain-containing RING finger protein 3) (SH3 multiple domains protein 4)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use