


Product Name: Recombinant Human TRIM38 Protein
Catalog #: PEH9614



Summary

Name	TRI38_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00635
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CGHSYCHLCI TDFFKNPSQK QLRQETFCCP QCRAPFHMDL LRPNKQLGSL IEALKETDQE MSCEEHGEQF HLFCEDEGQL ICWRCERAPQ HKGHTTALVE DVCQGYKEKL QKAVTKLKQL EDRCTEQKLS TAMRITKWKE KVQIQRQKIR SDFKNLQCFL HEEKSYLWR LEKEEQQTLS RLRDYEAGLG LKSNEKSHI LELEEKQGS AQKLLQNVND TLSRSWAVKL

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM38 (EC 2.3.2.27) (RING finger protein 15) (RING-type E3 ubiquitin transferase TRIM38) (Tripartite motif-containing protein 38) (Zinc finger protein RoRet)
--------------------------	--

Product Name: Recombinant Human TRIM38 Protein
Catalog #: PEH9614



Note

For Research Use Only , Not for Diagnostic Use