


Product Name: Recombinant Human TRIM7 Protein
Catalog #: PEH9615



Summary

Name	TRIM7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9C029
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQHGEPEFKLY CQDDGRAICV VCDRAREHRE HAVLPLDEAV QEAKELLESRLRVLKKELED CEVFRSTEKK ESKELLQMA AEQEKVGAEF QALRAFLVEQ EGRLLGRLEE LSREVAQKQN ENLAQLGVEI TQLSKLSSQI QETAQKPDLD

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM7 (EC 2.3.2.27) (Glycogenin-interacting protein) (RING finger protein 90) (Tripartite motif-containing protein 7)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use