


Product Name: Recombinant Human ELP4 Protein
Catalog #: PEH9618



Summary

Name	ELP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EB1
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YKDYHGLIHI RQIPRLNNLI CDESDVKDLA FKLKRKLFTI ERLHLPPDLS DTVSRSSKMD LAESAKRLGP GCGMMAGGKK

image



Background

Alternative Names Elongator complex protein 4 (hELP4) (PAX6 neighbor gene protein)

Note

For Research Use Only , Not for Diagnostic Use