


Product Name: Recombinant Human P3H4 Protein
Catalog #: PEH9619



Summary

Name	SC65_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92791
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CKVDCEANLT PNVGGYFVDK FVATMYHYLQ FAYYKLNDVR QAARSAASYM LFDPKDSVMQ QNLVYYRFHR ARWGLEEDF QPREEAMLYH NQTAEIRELL EFTHMYLQSD DEMELEETEP PLEPEDALSD AEFEGEGDYE EGMADWWQE PDAKGDEAEA

image



Background

Alternative Names	Endoplasmic reticulum protein SC65 (Leprecan-like protein 4) (Nucleolar autoantigen No55) (Prolyl 3-hydroxylase family member 4) (Synaptonemal complex protein SC65)
--------------------------	--

Note

Product Name: Recombinant Human P3H4 Protein
Catalog #: PEH9619



For Research Use Only , Not for Diagnostic Use