


Product Name: Recombinant Human EDC3 Protein
Catalog #: PEH9620



Summary

Name	EDC3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96F86
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSAPQNIPKR TDVKSQDVAV SPQQQCCKSKS YVDRHMESLS QSKSFRRRHNSWSSSRHPN QATPKKSGLK NGQMKNKDDE

image



Background

Alternative Names	Enhancer of mRNA-decapping protein 3 (LSM16 homolog) (YjeF N-terminal domain-containing protein 2) (YjeF_N2) (hYjeF_N2) (YjeF domain-containing protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use