


Product Name: Recombinant Human EIF3E Protein
Catalog #: PEH9622



Summary

Name	EIF3E_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P60228
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AHFLDRHLVF PLLEFLSVKE IYNEKELLQG KDLLSDTNM VDFAMDVYKN LYSDDIPHAL REKRTTVVAQ LKQLQAETEP IVKMFEDPET TRQMQSTRDG RMLFDYLADK HGFRQEYLDL LYRYAKFQYE CGNYSGAAEY LYFFRVLVPA TDRNALSSLW GKLASEILMQ NWDAAMEDLT RLKETIDNNS VSSPLQSLQQ

image



Background

Alternative Names	Eukaryotic translation initiation factor 3 subunit E (eIF3e) (Eukaryotic translation initiation factor 3 subunit 6) (Viral integration site protein INT-6 homolog) (eIF-3 p48)
--------------------------	--

Note

Product Name: Recombinant Human EIF3E Protein
Catalog #: PEH9622



For Research Use Only , Not for Diagnostic Use