


Product Name: Recombinant Human EIF3F Protein
Catalog #: PEH9623



Summary

Name	EIF3F_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O00303
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAPASSSDPA AAAAATAAPG QTPASQAQA QTPAPALPGP ALGPFPGGR VVRHPVILA SIVDSYERRN EGAARVIGTL LGTVDKHSVE VTNCFSVPHN ESEDEVAVDM EFAKNMYELH KKVSPNELIL GWYATGHDIT EHSVLIHEY SREAPNPIHL TVDTSLQNGR MSAKAYVSTL MGVPGRMTMGV MFTPLTVKYA YYDTERIGVD LIMKTCFSPN RVIGLSSDLQ QVGGASARIQ DALSTVLQYA EDVLSGKVSA DNTVGRFLMS LVNQVPKIVP DDFETMLNSN INDLLMVTYL ANLTQSQIAL

image



Background

Alternative Names	Eukaryotic translation initiation factor 3 subunit F (eIF3f) (Deubiquitinating enzyme eIF3f) (EC 3.4.19.12) (Eukaryotic translation initiation factor 3
--------------------------	---

Product Name: Recombinant Human EIF3F Protein
Catalog #: PEH9623



subunit 5) (eIF-3-epsilon) (eIF3 p47)

Note

For Research Use Only , Not for Diagnostic Use