


Product Name: Recombinant Human EIF4E1B Protein
Catalog #: PEH9625



Summary

Name	I4E1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NMX2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGIREWEEEE KEEAAERTP TGEKSPNSPR TLLSLRGKAR TGGPMEVKLE LHPLQNRWAL WFFKNDRSRA WQDNLHLVTK

image



Background

Alternative Names Eukaryotic translation initiation factor 4E type 1B

Note

For Research Use Only , Not for Diagnostic Use