


Product Name: Recombinant Human EXOC4 Protein
Catalog #: PEH9628



Summary

Name	EXOC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96A65
Host	E.coli
Species	Human
Molecular Mass	75 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKKPQRPKNS LKFESSHA ISMSAYLREQ RRELYSRSGE LQGGPDDNLI EGGGTKFVCK PGARNITVIF HPLLRFIQEI EHALGLGPAK QCPLREFLTV YIKNIFLNQV LAEINKEIEG VTKTSDPLKI LANADTMKVL GVQRPLLQST IIVEKTVQDL LNLMDLSAY SDQFLNMVCV KLQEYKDTCT AAYRGIVQSE EKLVISASWA KDDDISRLK SLPNWMNMAQ PKQLRPKREE EEDFIRAAFG KESEVLIGNL GDKLIPPQDI LRDVSDLKAL

image



Background

Alternative Names	Exocyst complex component 4 (Exocyst complex component Sec8)
--------------------------	--

Product Name: Recombinant Human EXOC4 Protein
Catalog #: PEH9628



Note

For Research Use Only , Not for Diagnostic Use