


Product Name: Recombinant Human EXOC4 Protein
Catalog #: PEH9629



Summary

Name	EXOC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q96A65
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEDFIRAAFG KESEVLIGNL GDKLIPPQDI LRDVSDLKAL ANMHESLEWL ASRTKSAFSN LSTSQMLSPA QDSHTNTDLP PVSEQIMQTL SELAKSFQDM ADRCLLVHL EVRVHCFHYL IPLAKEGNYA IVANVESMDY DPLVVKLNKD ISAIEEAMSA SLQQHKFQYI FEGLGHLISC ILINGAQYFR

image



Background

Alternative Names Exocyst complex component 4 (Exocyst complex component Sec8)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EXOC4 Protein
Catalog #: PEH9629

