


Product Name: Recombinant Human EXOSC3 Protein
Catalog #: PEH9630



Summary

Name	EXOS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQT5
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLAGSRARAA RTVLGQVVLP GEELLPEQE DAEGPGGAVE RPLSLNARAC SRVRVVC GPG LRRCGDRLLV TKCGRLRHKE PGSGSGGGVY WVDSQQKRYV PVKGDHVG I VTAKSGDIFK VDVGGSEPAS LSYLSFEGAT KRNRPNVQVG DLIYGQFVVA NKDMEPEMVC IDSCGRANGM GVIGQDGLLF KVTGLIRKL LAPDCEIIQE VGKLYPLEIV FGMNGRIWVK AKTIQQTLL ANILEACEHM

image



Background

Alternative Names	Exosome complex component RRP40 (Exosome component 3) (Ribosomal RNA-processing protein 40) (p10)
--------------------------	---

Note

Product Name: Recombinant Human EXOSC3 Protein
Catalog #: PEH9630



For Research Use Only , Not for Diagnostic Use