


**Product Name: Recombinant Human TBRG4 Protein**  
**Catalog #: PEH9634**



## Summary

<b>Name</b>	FAKD4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q969Z0
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	32 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	YHLVQKPFSL TKDVLLDVAY AYGKLSFHQT QVSQRLATDL LSLMPSLTSG EVAHCAKSFA LLKWLSLPLF EAFAQHVLNR AQDITLPHLC SVLLAFARLN FHPDQEDQFF SLVHEKLGSE LPGLEPALQV DLVWALCVLQ QAREAELQAV LHPEFHIQFL GGKSQKDQNT FQKLLHINAT ALLEYPEYSG PLLPASAVAP GPSALDRKVT PLQKELQETL KGLLSADKG SLEVATQYGW VLDAEVL LDS DGEFLPVRDF VAPHLAQPTG

## image



## Background

<b>Alternative Names</b>	FAST kinase domain-containing protein 4 (Cell cycle progression restoration protein 2) (Cell cycle progression protein 2) (Protein TBRG4) (Transforming growth factor beta regulator 4)
--------------------------	---

**Product Name: Recombinant Human TBRG4 Protein**  
**Catalog #: PEH9634**

---



**Note**

For Research Use Only , Not for Diagnostic Use